

Same DNA. More Solutions.

Valid from July 1, 2024 - September 30, 2024

FROM CANADA'S EXCLUSIVE AMANN GIRRBACH DISTRIBUTOR

MOTION 2 DRY & HYBRID PACKAGES

MOTION 2 DRY (Regular Price \$54,999) Now \$32,999



Item Code	Description	Retail Price	Quantity	Discount	Total Price
996902	CAM Software Dongle- Match 2	\$ 199	1		\$ 199
996002	Performance Plan- Match 2	\$ 2,499	1	\$ 200	\$ 2,299
INSTRAIN	In-person installation and training	\$ 999	1		\$ 999
SHIP	Shipping	\$ 1,999	1	\$ 500	\$ 1,499
761776	Zolid Fx Multi-Layer 14mm 0/A1	\$ 189	1	\$ 189	\$ -
761780	Zolid Fx Multi-Layer 14mm A2/A3	\$ 189	1	\$ 189	\$ -
761788	Zolid Fx Multi-Layer 14mm 0/B1	\$ 189	1	\$ 189	\$ -
761800	Zolid Fx Multi-Layer 14mm C1/C2	\$ 189	1	\$ 189	\$ -
761808	Zolid Fx Multi-Layer 14mm D2/D3	\$ 189	1	\$ 189	\$ -
MARGAY100	Green State sprue polisher	\$ 16	1	\$ 16	\$ -

MOTION 2 HYBRID (INCLUDES COOLSTREAM & TANK, Regular Price \$64,999) Now \$42,999

FEATURES

- 6 Tool spaces with an automatic tool changer
- Suction cup for increased suction performance and reduction of water spray
- Coolant supply tank connection allows for long operating cycles
- 5-Axis milling unit allows for wet & dry processing
- Super high frequency spindle; precise, stable & quiet running
- Tool length measuring unit; including fracture check & calibration



TOTAL MOTION 2 DRY PACKAGE \$37,995*

TOTAL MOTION 2 HYBRID PACKAGE \$47,995*

*Plus applicable taxes

BIFIX QM

BUY 3 BIFIX QM 10G SYRINGES, GET 1 FREE!

TOP SELLER!



- For all materials (metal, composite inlays and ceramic, also zirconium dioxide)
- Bubble-free perfect mixing
- Direct application
- Excellent adhesion to dentine, enamel and to metals
- Special one-step Ceramic Bond for optimal adhesion to porcelain
- Radiopaque
- Available in three shades
- Available as an all-in-one set



Sizzlin' summer savings from

keystone
industries

KeyDenture Try-In

- Cost effective for try-in dentures
- Exceptional surface quality
- Easy to characterize
- Polishes to high luster
- Biocompatible medical device
- Free of MMA, BPA, and phthalates

Available in A1, B1, and BL4



**FREE
OPTIMIZE
ISOPROPYL
ALCOHOL**

**WITH A PURCHASE OF
KEYDENTURE TRY IN !**



15% OFF!
All Keystone Lathe burs,
Hatho Lathe polishing
discs and Mizzy Heatless
wheels!

Versacryl

Thermo-Soft Adjustable Reline Heat Cure



- Full and partial denture relines
- Added lingual retention on lower dentures
- Obturator cases (cancer patients)
- Post dam areas on full uppers
- Gasket retention partial dentures
- Lingual retention on partial dentures

**15% OFF VERSACRYL
PRODUCTS WHILE
SUPPLIES LAST!**

MAGNUM HYPERONE TITANINUM DISCS



New product!



- Innovative material, grade 23 titanium
- Highly biocompatible, Low specific weight
- Optimal modulus of compatibility for a more stable result
- Available in 10,15,16, & 18mm discs

THICKNESS

PRICE

10MM	\$119.00
15MM	\$199.00
16MM	\$219.00
18MM	\$259.00

MESA
PRODUCER OF SPECIAL ALLOYS SINCE 1975

FRESH DEALS AND GREAT PRODUCTS FROM

HarvestDental

SURFACTANT



An alcohol-free solution is gentle on investments, ensuring precise cavitation & extra super sharp margins

8 oz \$41.00 #7060
32 oz \$84.00 #7061



CERAFIRM

A glycol-free liquid for use with high-fusing and low fusing ceramics. Ceramic stays where it's put!

250ml \$91.00 #201-0000
10ml \$66.000 #201-1000
kit \$109.00 #201-2000

MASTER DIE SEALER



Seals & protects. With MS1, the precious die is well protected and ready for die treatment.

1 oz \$39.00 #7169
8 oz \$69.00 #7176

PAIRS WELL!



MASTER DIE SEPARATOR

Now, isolate Wax, Ceramic and Composite on the master die with just one liquid!

1 oz \$39.00 #7180
8 oz \$69.00 #7183

10% OFF

HARVEST PRODUCTS WHEN YOU PURCHASE TWO OR MORE AT A TIME !

VERSAGLUE



Versaglu bonds all the materials you use in the laboratory, including gypsum, metal, acrylic and porcelain!

\$24.00
#2025 PX4 Thin
#2026 PX40 Medium

PAIRS WELL!



ZBASE

Designed to block out voids and undercuts, while retaining the ability to be molded.

Z Base block out material 10g \$21.00 #2150
Z Base Accelerator Spray 2oz \$15.00 #2151
Z Base Accelerator Spray refill with spray bottle \$59.00 #2152

PAIRS WELL!

DURA PATTERN

Dura Pattern is a one component pattern resin. Use for detailed cast & press patterns, attachments units and implant jigs.



\$109.00
#3000

SURE FIT DIE SPACERS



Sure Fit accurate cement space simulation in a best-in-class, thixotropic and scratch resistant varnish

\$33.00

Available in:
Blue #7151
Red #7158
Gold #7159
Silver #4160



PEEK BIOSOLUTION

20% OFF!

- Ideal for a number of indications; Crowns & Bridges, denture bases, occlusal splints, transversal connectors & sublingual brackets.
- Suitable for permanent, definitive & surgically invasive treatment
- Anti allergenic- Metal and Monomer free
- Flexural strength - 170 MPa (ISO-20795-1)
- Elastic Modulus - 5,100 MPa (white and A2/B2), 4,500 MPa (Gum)
- PEEK and titanium dioxide



ITEM NO.	DESCRIPTION	REGULAR PRICE	NOW
1018962	Peek 16mm white	\$ 449.00	\$ 359.00
1018963	Peek 20mm white	\$ 499.00	\$ 399.00
1018964	Peek 24mm white	\$ 539.00	\$ 431.00
1018967	Peek A2/B2 16mm	\$ 449.00	\$ 359.00
1018968	Peek A2/B2 20mm	\$ 499.00	\$ 399.00
1018969	Peek A2/B2 24mm	\$ 539.00	\$ 431.00
1019230	Peek Gum 16mm	\$ 449.00	\$ 359.00
1019231	Peek Gum 20mm	\$ 499.00	\$ 399.00
1019232	Peek Gum 24mm	\$ 539.00	\$ 431.00

AMANNGIRRBACH CERAMILL A-TEMP

- Highest esthetics due to perfectly matched VITA A-D shades for the Zolid DNA Generation
- Acrylic for long-term temporary restorations with a wearing period of up to three years
- Temporary anterior and posterior crowns

15% OFF

A1 CERAMILL A-TEMP A1 98X16

A3,5 CERAMILL A-TEMP A3.5 98X20

B2 CERAMILL A-TEMP B2 98X16

C2 CERAMILL A-TEMP C2 98X14

A2 CERAMILL A-TEMP A2 98X14

A3,5 CERAMILL A-TEMP A3,5 98X16

B2 CERAMILL A-TEMP B2 98X20



*ON SELECT PUCKS. ONLY WHILE SUPPLIES LAST

DETA SEAL HYDROFLOW SOFT PUTTY FROM DETAX

#02880

500ml jar of base
500ml jar of catalyst
2 measuring scoops



10% OFF



DETA SEAL PUTTY FROM DETA X

#02727

250ml jar of base
250ml jar of catalyst
2 measuring scoops

SAVE \$79 PER BOTTLE ON FREEPRINT SPLINT UV RESIN!



#03709

FLEXISTONE PLUS

An insulating & modelling material.
Great for Ivoclar or Ivobase system. Allows
for a protective barrier of acrylics & and
acrylic teeth against plaster.



#02383 2 X 160 ML, 1 MIXING BOWL, BASE + CATALYST

REGULAR PRICE \$139, NOW \$109. SAVE \$30

#02329 2 X 1000 ML, 4 MIXING BOWLS, BASE + CATALYST

REGULAR PRICE \$649, NOW \$519. SAVE \$130

Detaseal Hydroflow VPS



Regular Set

Working time: 2:00
Setting Time: 2:00

Fast Set

Working time: 1:00
Setting Time: 1:30

Mono Fast #02754 \$69.00


XLite, Fast #02744 \$69.00

XLite, Regular #02741 \$69.00

XLite, Regular #02741 \$69.00

XLite, Fast #02766 \$69.00

DETAX << PREMIUM BITE REGISTRATIONS



Glassbite

Working time 0:30

Setting time 1:00

#02574

\$39.00

Milkbite

Working time 0:30

Setting time 0:50

#03380

\$69.00

Greenbite Apple

Working time 0:30

Setting time 0:50

#02641

\$72.00

P R E C I S I O N B I T E M A T E R I A L S

Toll Free: 855-604-0465 • abutmentdirect.com/dentology • [f](https://www.facebook.com/dentologybyabutmentdirect) [i](https://www.instagram.com/dentologybyabutmentdirect) @dentologybyabutmentdirect

Promotions available from July 1, 2024 to September 30, 2024.

All prices exclude applicable taxes. Free shipping on orders over \$500 before taxes.